

Anti-PA14 IgG ELISA Protocol

Day #1

1. Bacterial Lysis and Protein Assay

- a. Transfer 1ml of bacterial sub-culture into 1.5ml Eppendorf
- b. Spin for 5' at 12,000 rpm
- c. Wash with 1 ml PBS; spin for 5' at 12,000 rpm
- d. Freeze and thaw the bacterial pellet ~3 times using dry ice
- e. Re-suspend the pellet in 100ul PBS and freeze thaw ~3 times
- f. Add another 900ul PBS, re-suspend well and measure protein concentration using Micro BCA Protein Assay Kit (dilute sample 1:100 or 1:10)

2. ELISA- Coating and Blocking

- a. Coat each well with 100ul PA14 lysate (10ug/ml, diluted in PBS) or PBS alone
- b. Incubate for 2hr at 37C
- c. Wash 3X with PBS + 0.05% Tween
- d. Add 300ul PBS+5% Normal Goat Serum (NGS) and incubate o/n at RT

Day #2

3. ELISA Development

- a. Wash coated ELISA plate 3X with PBS + 0.05% Tween
- b. Add 100ul patient serum diluted 1:100 in PBS+1% NGS
- c. Incubate for 1hr at RT
- d. Wash 3X with PBS + 0.05% Tween
- e. Add 100ul goat-anti-human HRP-labeled 2' ab (1:50,000 dil. in PBS+1% NGS) and incubate for 1hr at RT; take out Substrate Reagents and Stop soln. to warm up at RT, keep in the dark.
- f. Wash 3X with PBS + 0.05% Tween
- g. Add 100ul substrate (R&D, Substrate Reagent A: Substrate Reagent B 1:1)
- h. Incubate at RT in the dark for ~15min
- i. Add 50ul stop solution and read absorbance at 450nm and 570nm (for correction).